



Recharge Pulse Wave

www.ReChargePulseWave.com

Larry Basch, D.C., CCSP, CCEP, ICSC LBASCHDC@gmail.com Ph: 858-945-8118
EIN # 26-1657557 CA Resale Permit # 22-4191872 DUNS# 133874680 Cage # 3GN7



MODEL 2 ReCharge PEMF units are the **ONLY** high power PEMF manufacturer to be safety tested having UL / CL certification, manufactured in an ISO 13485 Certified Facility, is conducting University research, and is made in the USA. And this is the only brand approved for sale to local, state and federal agencies, Department of Defense, US Navy / Marine Corps. High Intensity PEMF creates an electric field within a magnetic coil and requires the use of a capacitor to store energy, which takes some time to build up. It will then discharge that stored energy in a microsecond. The frequency of discharge is related to the power setting machines. If you increase the power, the frequency is decreased due to the capacitor limits. Other low power units (often called "Hummers") use a high frequency oscillator with extremely low power output. The power output is so low; patients do not feel these fast oscillations passing through the body. This results in very long treatment times with little noticeable benefits. With our High Power PEMF, you feel it right away, so there is no doubt its working, you feel the pulse immediately deep within the body. This allows for the shortest treatment times & deeper tissue penetration.

ReCharge PEMF is the **ONLY Brand of PEMF** using Dual Coils delivering **Dual Pulsing** offering the proven benefits of both Low Power PEMF & High-Power PEMF running simultaneously so a second 10Hz baseline low power pulse is always running in addition to the higher intensity PEMF which you can select the power setting. So you get two different machines with two different pulses running simultaneously or (**DUAL PULSE FREQUENCY**) in one medical grade machine.

* (We offer three different Dual Pulse High Power units. This is the medium priced & medium power output of our high-power units, meaning slightly longer treatment times and slightly less depth of penetration as the highest power units.)

ReCharge Model #2 Dual Medium Mini Tote or Model #2 Dual Medium Clinical Desktop



Medium - High Power



3-year Warranty



ReCharge Model 2.0 BASIC PACKAGE Select 4 Attachments

ReCharge Model 2.1 COMPLETE CLINICAL PACKAGE Includes All 7 Attachments

Tote: Weight 9 lbs (12.5"W x 10.25D x 6"H)

Desktop: Weight 6.3 lbs (10"W x 10" D x 4"H)

Model 6 or 7 uses 2 SQUARE WAVES SIMULTANEOUSLY to carry the magnetic field directly to the cell at up to 5000 gauss. Instead of resonating just the receptors on the cell membrane, this Model delivers harmonic impulses that resonates the entire cell membrane. This is a high-intensity impulse system and is suitable for use both in-home or in a clinical setting. This unit a minimum gauss of 500 with a max of 5000, but fully adjustable intensity.

ReCharge Attachments

Single Loop(10") / Small Butterfly Loop (7")
Paddle (6") / Pelvic Cushion / Brain Helmet
Small Mat (12"x22") / Large Body Mat (24" x 72")
Versatile Use: Home - Clinic - Travel
5-10 minute Timer

Low Mode Pulses/Sec: 10.
High Mode Pulses/Sec: 1 + 10 (Low).
Power High Mode Intensity: estimate **5220 Gauss** =.5 Tesla.
Low Intensity: estimate 400 Gauss
Voltage Regulator: Auto adapts to 120-230 V,
50-60 Hz AC, worldwide.

MODEL 2 Tote or DESKTOP

Equipment Purchase / Finance

More Detailed PEMF info at: www.ReChargePulseWave.com LBASCHDC@gmail.com Ph: 858-945-8118

Stearns Bank	Ellen Boquist	800-247-1922 Ext. 3199	cell: 320-248-7116	EllenB@stearnsbank.com
NCMIC	Scott Fuller	877-770-7244, ext. 4685	cell 515-313-4685	sfuller@ncmic.com ncmic.com